

 amplify your imaging value™



# iConnect® Gateway V. 1.0

## DICOM CONFORMANCE STATEMENT

Merge Healthcare  
900 Walnut Ridge Drive  
Hartland, WI 53029  
USA

**Copyright 2016 Merge Healthcare, an IBM Company**  
**Unauthorized use, reproduction, or disclosure is prohibited.**

Use, duplication, or disclosure by licensee is subject to restrictions as set forth in the Merge Healthcare contract. Information contained in this document is subject to change without notice and does not represent a commitment on the part of Merge Healthcare. The software and processes described in this document are furnished under the contract or nondisclosure agreement. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording, for any purpose other than the purchaser's personal use without the written permission of Merge Healthcare.

This document contains proprietary information of Merge Healthcare. The contents are exempt from disclosure to the public under the Freedom of Information Act 35, U.S.C. 552 (6)(4) and unlawful disclosure thereof is a violation of the Trade Secrets Act, 18 U.S.C. 1905. Public disclosure of any information contained herein shall not be made without prior written permission of Merge Healthcare.

iConnect® and Merge Healthcare® are registered trademarks of Merge Healthcare, Inc.

The Merge Healthcare logo is a trademark of Merge Healthcare.

All other names and products are trademarks or registered trademarks of their respective companies.

### **U.S. Government Restricted Rights**

This product is a "Commercial Item" offered with "Restricted Rights." The Government's rights to use, modify, reproduce, release, perform, display or disclose this documentation are subject to the restrictions set forth in Federal Acquisition Regulation ("FAR") 12.211 and 12.212 for civilian agencies and in DFARS 227.7202-3 for military agencies. Contractor is Merge Healthcare.

The information contained herein is confidential and is the sole property of Merge Healthcare. None of the information may be reproduced, copied, and/or re-distributed electronically, manually, or otherwise to any person without advance written permission from Merge Healthcare, except as may be expressly authorized by contracts or other agreements between Merge Healthcare and its customer. Merge Healthcare® and AMICAS are registered trademarks and VERICIS, RadSuite, Abraxas, and iConnect are also trademarks of Merge Healthcare. Other trademarks or service marks referred to herein are the property of their respective owners.

For application support or to report issues with user documentation, contact Customer Support:

- In North America: call toll-free 1-877-741-5369
- Outside of North America: +31.40.299.0773
- Email: support@merge.com

| Part     | Date          | Revision | Description  |
|----------|---------------|----------|--|
| CLD-2928 | November 2014 | 1.0      | Released under iConnect Cloud Project.               |
| IG-5     | October 2015  | 2.0      | Released under iConnect Gateway Project. No changes. |
|          | January 2016  | 3.0      | Update to release date. No changes.                  |

# Contents

---

|   |           |
|---|-----------|
| <b>Contents</b>                                     | <b>i</b>  |
| <b>Introduction</b>                                 | <b>4</b>  |
| Purpose .....                                       | 4         |
| Related Documents .....                             | 4         |
| <b>Implementation Model</b> .....                   | <b>5</b>  |
| iConnect Gateway .....                              | 5         |
| Application Data Flow Diagram .....                 | 5         |
| Functional Definition of AEs .....                  | 6         |
| <b>AE Specifications</b> .....                      | <b>7</b>  |
| iConnect Gateway DICOM Specifications .....         | 7         |
| Association Establishment Policies .....            | 10        |
| Association Initiation by Real-World Activity ..... | 10        |
| Association Acceptance Policy .....                 | 15        |
| <b>Configuration</b>                                | <b>32</b> |
| AE Title/Presentation Address Mapping .....         | 32        |
| Configurable Parameters .....                       | 32        |

---

---

# Introduction

---

## Purpose

This document specifies the conformance of Merge Healthcare's iConnect Gateway Release V.1.1.0 to the DICOM 3.0 standard. This document covers conformance for the following:

- Merge iConnect Gateway – a DICOM to HTTPS transition server used for querying and transferring DICOM images between Other DICOM devices and the Merge iConnect Cloud Archive.

---

## Related Documents


NEMA, the DICOM Standard: Parts 1 - 20.

# Implementation Model

## iConnect Gateway

The iConnect Gateway allows for the storage and retrieval of DICOM objects from the iConnect Cloud Archive. Additionally, it allows for the querying of information about the DICOM objects that have been stored to it the iConnect Cloud Archive.

### Application Data Flow Diagram



iConnect Gateway DICOM AE Implementation Model

---

## Functional Definition of AEs

---

---

### DICOM AE

---

The iConnect Gateway waits for another application to connect at the presentation address configured for its Application Entity Title. When another application connects, the archive expects it to be a DICOM application. The archive will accept associations with Presentation Contexts for the following SOP Classes:

- Storage Service Class
- Query/Retrieve Service Class

The archive will receive objects/requests. When it receives an object, it stores that object in dicom Part 10 format as a file on a file system and registers that object in the database.

When a request is received it is processed accordingly. If the request is a retrieve request, the archive will retrieve the object from Cloud Archive in which it is stored. Once retrieved, it will be sent to the desired destination via an association. . If the request is a query, the archive will simply perform the query to the iConnect Cloud Archive and return the results.

# AE Specifications

## iConnect Gateway DICOM Specifications

This Application Entity provides Standard Conformance to the following DICOM V3.0 SOP Classes as an SCU and/or SCP (as indicated in the table):

**NOTE:** New SOP classes may be supported through simple configuration changes of the product.

**Table 1. SOP Classes**

| SOP Class Name  | SOP Class UID                 | SCU/SCP |
|---|-------------------------------|---------|
| <b>Storage</b>  |                               |         |
| 12-lead ECG Waveform Storage                          | 1.2.840.10008.5.1.4.1.1.9.1.1 | Y / Y   |
| Ambulatory ECG Waveform Storage                       | 1.2.840.10008.5.1.4.1.1.9.1.3 | Y / Y   |
| Basic Text SR Storage                                 | 1.2.840.10008.5.1.4.1.1.88.11 | Y / Y   |
| Basic Voice Audio Waveform Storage                    | 1.2.840.10008.5.1.4.1.1.9.4.1 | Y / Y   |
| Cardiac Electrophysiology Waveform Storage            | 1.2.840.10008.5.1.4.1.1.9.3.1 | Y / Y   |
| Comprehensive SR Storage                              | 1.2.840.10008.5.1.4.1.1.88.33 | Y / Y   |
| Computed Radiography Image Storage                    | 1.2.840.10008.5.1.4.1.1.1     | Y / Y   |
| CT Image Storage                                      | 1.2.840.10008.5.1.4.1.1.2     | Y / Y   |
| Enhanced CT Image Storage                             | 1.2.840.10008.5.1.4.1.1.2.1   | Y / Y   |
| Digital X-Ray Image Storage (Presentation)            | 1.2.840.10008.5.1.4.1.1.1.1   | Y / Y   |
| Digital X-Ray Image Storage (Processing)              | 1.2.840.10008.5.1.4.1.1.1.1.1 | Y / Y   |
| Digital Mammography Image Storage (Presentation)      | 1.2.840.10008.5.1.4.1.1.1.2   | Y / Y   |
| Digital Mammography Image Storage (Processing)        | 1.2.840.10008.5.1.4.1.1.1.2.1 | Y / Y   |
| Digital Intra-oral X-Ray Image Storage (Presentation) | 1.2.840.10008.5.1.4.1.1.1.3   | Y / Y   |
| Digital Intra-oral X-Ray Image Storage (Processing)   | 1.2.840.10008.5.1.4.1.1.1.3.1 | Y / Y   |
| Encapsulated PDF Storage SOP Class                    | 1.2.840.10008.5.1.4.1.1.104.1 | Y / Y   |
| Enhanced SR Storage                                   | 1.2.840.10008.5.1.4.1.1.88.22 | Y / Y   |

|  |                                  |       |
|--|----------------------------------|-------|
| General ECG Waveform Storage                   | 1.2.840.10008.5.1.4.1.1.9.1.2    | Y / Y |
| Grayscale Softcopy Presentation State Storage  | 1.2.840.10008.5.1.4.1.1.11.1     | Y / Y |
| Hardcopy Color Image Storage                   | 1.2.840.10008.5.1.1.30           | Y / Y |
| Hardcopy Grayscale Image Storage               | 1.2.840.10008.5.1.1.29           | Y / Y |
| Hemodynamic Waveform Storage                   | 1.2.840.10008.5.1.4.1.1.9.2.1    | Y / Y |
| Key Object Selection Document                  | 1.2.840.10008.5.1.4.1.1.88.59    | Y / Y |
| Mammography CAD SR                             | 1.2.840.10008.5.1.4.1.1.88.50    | Y / Y |
| Multi-frame Single Bit Secondary Capture Image | 1.2.840.10008.5.1.4.1.1.7.1      | Y / Y |
| Multi-frame Grayscale Byte Secondary Capture   | 1.2.840.10008.5.1.4.1.1.7.2      | Y / Y |
| Multi-frame Grayscale Word Secondary Capture   | 1.2.840.10008.5.1.4.1.1.7.3      | Y / Y |
| Multi-frame True Color Secondary Capture Image | 1.2.840.10008.5.1.4.1.1.7.4      | Y / Y |
| MR Image Storage                               | 1.2.840.10008.5.1.4.1.1.4        | Y / Y |
| Enhanced MR Image Storage                      | 1.2.840.10008.5.1.4.1.1.4.1      | Y / Y |
| Nuclear Medicine Image Storage                 | 1.2.840.10008.5.1.4.1.1.20       | Y / Y |
| Nuclear Medicine Image Storage (Retired)       | 1.2.840.10008.5.1.4.1.1.5        | Y / Y |
| Ophthalmic Photography 8 Bit Image Storage     | 1.2.840.10008.5.1.4.1.1.77.1.5.1 | Y / Y |
| Ophthalmic Photography 16 Bit Image Storage    | 1.2.840.10008.5.1.4.1.1.77.1.5.2 | Y / Y |
| Ophthalmic Tomography Image Storage            | 1.2.840.10008.5.1.4.1.1.77.1.5.4 | Y / Y |
| Positron Emission Tomography Image Storage     | 1.2.840.10008.5.1.4.1.1.128      | Y / Y |
| RT Beams Treatment Record Storage              | 1.2.840.10008.5.1.4.1.1.481.4    | Y / Y |
| RT Brachy Treatment Record Storage             | 1.2.840.10008.5.1.4.1.1.481.6    | Y / Y |
| RT Dose Storage                                | 1.2.840.10008.5.1.4.1.1.481.2    | Y / Y |
| RT Image Storage                               | 1.2.840.10008.5.1.4.1.1.481.1    | Y / Y |
| RT Plan Storage                                | 1.2.840.10008.5.1.4.1.1.481.5    | Y / Y |
| RT Structure Set Storage                       | 1.2.840.10008.5.1.4.1.1.481.3    | Y / Y |
| RT Treatment Summary Record Storage            | 1.2.840.10008.5.1.4.1.1.481.7    | Y / Y |
| Secondary Capture Image Storage                | 1.2.840.10008.5.1.4.1.1.7        | Y / Y |
| Stand-alone Curve Storage                      | 1.2.840.10008.5.1.4.1.1.9        | Y / Y |



|  |                                |       |
|--|--------------------------------|-------|
| Stand-alone Modality LUT Storage                           | 1.2.840.10008.5.1.4.1.1.10     | Y / Y |
| Stand-alone Overlay Storage                                | 1.2.840.10008.5.1.4.1.1.8      | Y / Y |
| Stand-alone VOI LUT Storage                                | 1.2.840.10008.5.1.4.1.1.11     | Y / Y |
| Standalone PET Curve Storage                               | 1.2.840.10008.5.1.4.1.1.129    | Y / Y |
| Stored Print Storage                                       | 1.2.840.10008.5.1.1.27         | Y / Y |
| Ultrasound Image Storage                                   | 1.2.840.10008.5.1.4.1.1.6.1    | Y / Y |
| Ultrasound Image Storage (Retired)                         | 1.2.840.10008.5.1.4.1.1.6      | Y / Y |
| Ultrasound Multi-frame Image Storage                       | 1.2.840.10008.5.1.4.1.1.3.1    | Y / Y |
| Ultrasound Multi-frame Image Storage (Retired)             | 1.2.840.10008.5.1.4.1.1.3      | Y / Y |
| VL Endoscopic Image Storage                                | 1.2.840.10008.5.1.4.1.1.77.1.1 | Y / Y |
| VL Microscopic Image Storage                               | 1.2.840.10008.5.1.4.1.1.77.1.2 | Y / Y |
| VL Slide-Coordinates Microscopic Image Storage             | 1.2.840.10008.5.1.4.1.1.77.1.3 | Y / Y |
| VL Photographic Image Storage                              | 1.2.840.10008.5.1.4.1.1.77.1.4 | Y / Y |
| X-Ray Angiographic Bi-Plane Image Storage (retired)        | 1.2.840.10008.5.1.4.1.1.12.3   | Y / Y |
| X-Ray Angiographic Image Storage                           | 1.2.840.10008.5.1.4.1.1.12.1   | Y / Y |
| X-Ray Radiofluoroscopic Image Storage                      | 1.2.840.10008.5.1.4.1.1.12.2   | Y / Y |
| <b>Query/Retrieve</b>                                      |                                |       |
| Patient Root Query/Retrieve Information Model - FIND       | 1.2.840.10008.5.1.4.1.2.1.1    | N / Y |
| Patient Root Query/Retrieve Information Model - MOVE       | 1.2.840.10008.5.1.4.1.2.1.2    | N / Y |
| Patient Root Query/Retrieve Information Model - GET        | 1.2.840.10008.5.1.4.1.2.1.3    | N / Y |
| Study Root Query/Retrieve Information Model - FIND         | 1.2.840.10008.5.1.4.1.2.2.1    | N / Y |
| Study Root Query/Retrieve Information Model - MOVE         | 1.2.840.10008.5.1.4.1.2.2.2    | N / Y |
| Study Root Query/Retrieve Information Model - GET          | 1.2.840.10008.5.1.4.1.2.2.3    | N / Y |
| Patient/Study Only Query/Retrieve Information Model - FIND | 1.2.840.10008.5.1.4.1.2.3.1    | N / Y |
| Patient/Study Only Query/Retrieve Information Model - MOVE | 1.2.840.10008.5.1.4.1.2.3.2    | N / Y |

|   |                             |       |
|---|-----------------------------|-------|
| Patient/Study Only Query/Retrieve Information Model - GET | 1.2.840.10008.5.1.4.1.2.3.3 | N / Y |
|---|-----------------------------|-------|

---

## Association Establishment Policies

---

### General

---

The DICOM interface of the iConnect Gateway will attempt to establish an association anytime a C-MOVE request is received from a remote application entity in order to store the requested data in the move. The archive will only attempt to establish associations in response to valid C-MOVE requests for images that are known to the iConnect Cloud Archive database.

The gateway will only attempt to C-Move to know destinations that must be configured before attempting a move.

The DICOM Component supports SCU/SCP Role Negotiation.

The DICOM Component supports configurable Maximum PDU sizes, both sent and received.

---

### Number of Associations

---

The iConnect Gateway will accept any number of simultaneous associations, the number of which can be configured. The maximum limit on the number of simultaneous associations is dependent on the number of open file descriptors allowed by the underlying operating system though the practical count will be limited by system resources and network bandwidth.

---

### Association Initiation by Real-World Activity

---

The iConnect Gateway will attempt to initiate one storage association for each C-MOVE request request that is received.

---

### Storage Association

---

- Completing a Move Request to a Remote System and
- Unsolicited Storage of DICOM Objects to the iConnect Cloud Archive

*Associated Real-World Activity*

The associated real-world activity is the receipt of a C-MOVE request from a remote application entity.

*Proposed Presentation Contexts*

The iConnect Gateway will propose a collection of presentation contexts, which will be obtained by applying the algorithm presented in the next section to the following Presentation Context Table:

**Table 2. Presentation Context Table**

| Abstract Syntax                                  |                                | Transfer Syntax             | Role | Extended Negotiation |
|--|--------------------------------|-----------------------------|------|----------------------|
| Name   | UID                            |                             |      |                      |
| 12-lead ECG Waveform Storage                     | 1.2.840.10008.5.1.4.1.1.9.1.1  | See <a href="#">Table 3</a> | SCU  | None                 |
| Ambulatory ECG Waveform Storage                  | 1.2.840.10008.5.1.4.1.1.9.1.3  | See <a href="#">Table 3</a> | SCU  | None                 |
| Basic Study Content Notification                 | 1.2.840.10008.1.9              | See <a href="#">Table 3</a> | SCU  | None                 |
| Basic Text SR Storage                            | 1.2.840.10008.5.1.4.1.1.88.11  | See <a href="#">Table 3</a> | SCU  | None                 |
| Basic Voice Audio Waveform Storage               | 1.2.840.10008.5.1.4.1.1.9.4.1  | See <a href="#">Table 3</a> | SCU  | None                 |
| Breast Tomosynthesis Image Storage               | 1.2.840.10008.5.1.4.1.1.13.1.3 | See <a href="#">Table 3</a> | SCU  | None                 |
| Cardiac Electrophysiology Waveform Storage       | 1.2.840.10008.5.1.4.1.1.9.3.1  | See <a href="#">Table 3</a> | SCU  | None                 |
| Comprehensive SR Storage                         | 1.2.840.10008.5.1.4.1.1.88.33  | See <a href="#">Table 3</a> | SCU  | None                 |
| Computed Radiography Image Storage               | 1.2.840.10008.5.1.4.1.1.1      | See <a href="#">Table 3</a> | SCU  | None                 |
| CT Image Storage                                 | 1.2.840.10008.5.1.4.1.1.2      | See <a href="#">Table 3</a> | SCU  | None                 |
| Enhanced CT Image Storage                        | 1.2.840.10008.5.1.4.1.1.2.1    | See <a href="#">Table 3</a> | SCU  | None                 |
| Digital X-Ray Image Storage (Presentation)       | 1.2.840.10008.5.1.4.1.1.1.1    | See <a href="#">Table 3</a> | SCU  | None                 |
| Digital X-Ray Image Storage (Processing)         | 1.2.840.10008.5.1.4.1.1.1.1.1  | See <a href="#">Table 3</a> | SCU  | None                 |
| Digital Mammography Image Storage (Presentation) | 1.2.840.10008.5.1.4.1.1.1.2    | See <a href="#">Table 3</a> | SCU  | None                 |

|  |                               |                             |     |      |
|--|-------------------------------|-----------------------------|-----|------|
| Digital Mammography Image Storage (Processing)             | 1.2.840.10008.5.1.4.1.1.1.2.1 | See <a href="#">Table 3</a> | SCU | None |
| Digital Intra-oral X-Ray Image Storage (Presentation)      | 1.2.840.10008.5.1.4.1.1.1.3   | See <a href="#">Table 3</a> | SCU | None |
| Digital Intra-oral X-Ray Image Storage (Processing)        | 1.2.840.10008.5.1.4.1.1.1.3.1 | See <a href="#">Table 3</a> | SCU | None |
| Encapsulated PDF Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.104.1 | See <a href="#">Table 3</a> | SCU | None |
| Enhanced SR Storage  | 1.2.840.10008.5.1.4.1.1.88.22 | See <a href="#">Table 3</a> | SCU | None |
| General ECG Waveform Storage                               | 1.2.840.10008.5.1.4.1.1.9.1.2 | See <a href="#">Table 3</a> | SCU | None |
| Grayscale Softcopy Presentation State Storage              | 1.2.840.10008.5.1.4.1.1.11.1  | See <a href="#">Table 3</a> | SCU | None |
| Hardcopy Color Image Storage                               | 1.2.840.10008.5.1.1.30        | See <a href="#">Table 3</a> | SCU | None |
| Hardcopy Grayscale Image Storage                           | 1.2.840.10008.5.1.1.29        | See <a href="#">Table 3</a> | SCU | None |
| Hemodynamic Waveform Storage                               | 1.2.840.10008.5.1.4.1.1.9.2.1 | See <a href="#">Table 3</a> | SCU | None |
| Key Object Selection Document                              | 1.2.840.10008.5.1.4.1.1.88.59 | See <a href="#">Table 3</a> | SCU | None |
| Mammography CAD SR   | 1.2.840.10008.5.1.4.1.1.88.50 | See <a href="#">Table 3</a> | SCU | None |
| Multi-frame Single Bit Secondary Capture Image Storage     | 1.2.840.10008.5.1.4.1.1.7.1   | See <a href="#">Table 3</a> | SCU | None |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2   | See <a href="#">Table 3</a> | SCU | None |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3   | See <a href="#">Table 3</a> | SCU | None |
| Multi-frame True Color Secondary                           | 1.2.840.10008.5.1.4.1.1.7.4   | See <a href="#">Table 3</a> | SCU | None |
| MR Image Storage   | 1.2.840.100008.5.1.4.1.1.4    | See <a href="#">Table 3</a> | SCU | None |
| Enhanced MR Image Storage                                  | 1.2.840.100008.5.1.4.1.1.4.1  | See <a href="#">Table 3</a> | SCU | None |
| Nuclear Medicine Image Storage                             | 1.2.840.10008.5.1.4.1.1.20    | See <a href="#">Table 3</a> | SCU | None |

|   |                                  |                             |     |      |
|---|----------------------------------|-----------------------------|-----|------|
| Nuclear Medicine Image Storage              | 1.2.840.10008.5.1.4.1.1.5        | See <a href="#">Table 3</a> | SCU | None |
| Ophthalmic Photography 8 Bit Image Storage  | 1.2.840.10008.5.1.4.1.1.77.1.5.1 | See <a href="#">Table 3</a> | SCU | None |
| Ophthalmic Photography 16 Bit Image Storage | 1.2.840.10008.5.1.4.1.1.77.1.5.2 | See <a href="#">Table 3</a> | SCU | None |
| Ophthalmic Tomography Image Storage         | 1.2.840.10008.5.1.4.1.1.77.1.5.4 | See <a href="#">Table 3</a> | SCU | None |
| Positron Emission Tomography Image Storage  | 1.2.840.10008.5.1.4.1.1.128      | See <a href="#">Table 3</a> | SCU | None |
| RT Beams Treatment Record Storage           | 1.2.840.10008.5.1.4.1.1.481.4    | See <a href="#">Table 3</a> | SCU | None |
| RT Brachy Treatment Record Storage          | 1.2.840.10008.5.1.4.1.1.481.6    | See <a href="#">Table 3</a> | SCU | None |
| RT Dose Storage                             | 1.2.840.10008.5.1.4.1.1.481.2    | See <a href="#">Table 3</a> | SCU | None |
| RT Image Storage                            | 1.2.840.10008.5.1.4.1.1.481.1    | See <a href="#">Table 3</a> | SCU | None |
| RT Plan Storage                             | 1.2.840.10008.5.1.4.1.1.481.5    | See <a href="#">Table 3</a> | SCU | None |
| RT Structure Set Storage                    | 1.2.840.10008.5.1.4.1.1.481.3    | See <a href="#">Table 3</a> | SCU | None |
| RT Treatment Summary Record                 | 1.2.840.10008.5.1.4.1.1.481.7    | See <a href="#">Table 3</a> | SCU | None |
| Secondary Capture Image Storage             | 1.2.840.10008.5.1.4.1.1.7        | See <a href="#">Table 3</a> | SCU | None |
| Stand-alone Curve Storage                   | 1.2.840.10008.5.1.4.1.1.9        | See <a href="#">Table 3</a> | SCU | None |
| Stand-alone Modality LUT Storage            | 1.2.840.10008.5.1.4.1.1.10       | See <a href="#">Table 3</a> | SCU | None |
| Stand-alone Overlay Storage                 | 1.2.840.10008.5.1.4.1.1.8        | See <a href="#">Table 3</a> | SCU | None |
| Stand-alone VOI LUT Storage                 | 1.2.840.10008.5.1.4.1.1.11       | See <a href="#">Table 3</a> | SCU | None |
| Standalone PET Curve Storage                | 1.2.840.10008.5.1.4.1.1.129      | See <a href="#">Table 3</a> | SCU | None |
| Stored Print Storage                        | 1.2.840.10008.5.1.1.27           | See <a href="#">Table 3</a> | SCU | None |
| Ultrasound Image Storage                    | 1.2.840.10008.5.1.4.1.1.6.1      | See <a href="#">Table 3</a> | SCU | None |
| Ultrasound Image Storage (Retired)          | 1.2.840.10008.5.1.4.1.1.6        | See <a href="#">Table 3</a> | SCU | None |

|                                   |                                |                             |     |      |
|-----------------------------------|--------------------------------|-----------------------------|-----|------|
| Ultrasound Multi-frame Image      | 1.2.840.10008.5.1.4.1.1.3.1    | See <a href="#">Table 3</a> | SCU | None |
| Ultrasound Multi-frame Image      | 1.2.840.10008.5.1.4.1.1.3      | See <a href="#">Table 3</a> | SCU | None |
| VL Endoscopic Image Storage       | 1.2.840.10008.5.1.4.1.1.77.1.1 | See <a href="#">Table 3</a> | SCU | None |
| VL Microscopic Image Storage      | 1.2.840.10008.5.1.4.1.1.77.1.2 | See <a href="#">Table 3</a> | SCU | None |
| VL Slide-Coordinates Microscopic  | 1.2.840.10008.5.1.4.1.1.77.1.3 | See <a href="#">Table 3</a> | SCU | None |
| VL Photographic Image Storage     | 1.2.840.10008.5.1.4.1.1.77.1.4 | See <a href="#">Table 3</a> | SCU | None |
| X-Ray Angiographic Bi-Plane Image | 1.2.840.10008.5.1.4.1.1.12.3   | See <a href="#">Table 3</a> | SCU | None |
| X-Ray Angiographic Image Storage  | 1.2.840.10008.5.1.4.1.1.12.1   | See <a href="#">Table 3</a> | SCU | None |
| X-Ray Radiofluoroscopic Image     | 1.2.840.10008.5.1.4.1.1.12.2   | See <a href="#">Table 3</a> | SCU | None |

**Table 3. Transfer Syntax for Send to a Remote System**

| Name  | UID                    |
|---|------------------------|
| Little Endian Explicit VR                               | 1.2.840.10008.1.2.1    |
| Little Endian Implicit VR                               | 1.2.840.10008.1.2      |
| Big Endian Explicit VR                                  | 1.2.840.10008.1.2.2    |
| Deflated Little Endian Explicit VR                      | 1.2.840.10008.1.2.1.99 |
| Lossless JPEG Image Compression (baseline)              | 1.2.840.10008.1.2.4.70 |
| Lossy JPEG Image Compression (8-bit, coding Process 1)  | 1.2.840.10008.1.2.4.50 |
| Lossy JPEG Image Compression (12-bit, coding Process 4) | 1.2.840.10008.1.2.4.51 |
| JPEG Lossless, Nonhierarchical (Processes 14)           | 1.2.840.10008.1.2.4.57 |
| JPEG 2000 Image Compression (Lossless Only)             | 1.2.840.10008.1.2.4.90 |
| JPEG 2000 Image Compression                             | 1.2.840.10008.1.2.4.91 |
| Lossless RLE Image Compression                          | 1.2.840.10008.1.2.5    |

---

## Algorithm for Computing the Proposed Presentation Contexts for Outgoing Storage Associations

---

At the time of the association establishment, the set of DICOM objects to be sent out has been determined. Since both the transfer syntax in which each DICOM object has been stored and the SOP Class UID of each DICOM object are stored in the archives database, both of these tags are already known. The intention of the algorithm is to send out each object in the transfer syntax it has been initially stored in or a similar one. So, we understand that if the object has been initially received in a supported encapsulated transfer syntax, the algorithm would make sure that for the specific abstract syntax, a presentation context with that particular transfer syntax has been proposed. The encapsulated transfer syntax is preserved for efficiency, in order to save the processing time of the compression conversion. In the case of objects stored in a native transfer syntax, preserving the transfer syntax is not important, and any native transfer syntax should be acceptable.

---

## Association Acceptance Policy

When the iConnect Gateway accepts an association, it will allow the storage of objects, the retrieval of objects previously stored, the query for information about stored objects.

The iConnect Gateway is configurable to allow security restrictions on specified SOP classes, SCU/SCP roles and proposed presentation contexts. However, it accepts associations from any remote application entities.

---

## Receive Objects from a Remote System

---

### *Associated Real-World Activity*

When the iConnect Gateway receives a storage request (C-STORE) for an object in the iConnect Cloud Archive, that object is temporarily stored on disk on the underlying platform in DICOM Part 10 format and then upload into the iConnect Cloud Archive. If it is unable to store the object a failure response will be returned and the object will be marked as an error and DICOM object remained on disk..

The data set of the C-STORE command is stored with no loss of information.

### *Presentation Context Table*

The iConnect Gateway will accept the following Transfer Syntaxes for each presentation context in an association it receives:

**Table 4. Presentation Context Table**

| Abstract Syntax                                       |                                | Transfer Syntax             | Role | Extended Negotiation |
|---|--------------------------------|-----------------------------|------|----------------------|
| Name  | UID                            |                             |      |                      |
| 12-lead ECG Waveform Storage                          | 1.2.840.10008.5.1.4.1.1.9.1.1  | See <a href="#">Table 5</a> | SCP  | None                 |
| Ambulatory ECG Waveform Storage                       | 1.2.840.10008.5.1.4.1.1.9.1.3  | See <a href="#">Table 5</a> | SCP  | None                 |
| Basic Study Content Notification                      | 1.2.840.10008.1.9              | See <a href="#">Table 5</a> | SCP  | None                 |
| Basic Text SR Storage                                 | 1.2.840.10008.5.1.4.1.1.88.11  | See <a href="#">Table 5</a> | SCP  | None                 |
| Basic Voice Audio Waveform Storage                    | 1.2.840.10008.5.1.4.1.1.9.4.1  | See <a href="#">Table 5</a> | SCP  | None                 |
| Breast Tomosynthesis Image Storage                    | 1.2.840.10008.5.1.4.1.1.13.1.3 | See <a href="#">Table 5</a> | SCP  | None                 |
| Cardiac Electrophysiology Waveform Storage            | 1.2.840.10008.5.1.4.1.1.9.3.1  | See <a href="#">Table 5</a> | SCP  | None                 |
| Comprehensive SR Storage                              | 1.2.840.10008.5.1.4.1.1.88.33  | See <a href="#">Table 5</a> | SCP  | None                 |
| Computed Radiography Image Storage                    | 1.2.840.10008.5.1.4.1.1.1      | See <a href="#">Table 5</a> | SCP  | None                 |
| CT Image Storage                                      | 1.2.840.10008.5.1.4.1.1.2      | See <a href="#">Table 5</a> | SCP  | None                 |
| Enhanced CT Image Storage                             | 1.2.840.10008.5.1.4.1.1.2.1    | See <a href="#">Table 5</a> | SCP  | None                 |
| Digital X-Ray Image Storage (Presentation)            | 1.2.840.10008.5.1.4.1.1.1.1    | See <a href="#">Table 5</a> | SCP  | None                 |
| Digital X-Ray Image Storage (Processing)              | 1.2.840.10008.5.1.4.1.1.1.1.1  | See <a href="#">Table 5</a> | SCP  | None                 |
| Digital Mammography Image Storage (Presentation)      | 1.2.840.10008.5.1.4.1.1.1.2    | See <a href="#">Table 5</a> | SCP  | None                 |
| Digital Mammography Image Storage (Processing)        | 1.2.840.10008.5.1.4.1.1.1.2.1  | See <a href="#">Table 5</a> | SCP  | None                 |
| Digital Intra-oral X-Ray Image Storage (Presentation) | 1.2.840.10008.5.1.4.1.1.1.3    | See <a href="#">Table 5</a> | SCP  | None                 |
| Digital Intra-oral X-Ray Image Storage (Processing)   | 1.2.840.10008.5.1.4.1.1.1.3.1  | See <a href="#">Table 5</a> | SCP  | None                 |
| Encapsulated PDF Storage SOP Class                    | 1.2.840.10008.5.1.4.1.1.104.1  | See <a href="#">Table 5</a> | SCP  | None                 |
| Enhanced SR Storage                                   | 1.2.840.10008.5.1.4.1.1.88.22  | See <a href="#">Table 5</a> | SCP  | None                 |
| General ECG Waveform Storage                          | 1.2.840.10008.5.1.4.1.1.9.1.2  | See <a href="#">Table 5</a> | SCP  | None                 |
| Grayscale Softcopy Presentation State Storage         | 1.2.840.10008.5.1.4.1.1.11.1   | See <a href="#">Table 5</a> | SCP  | None                 |



|  |                                  |                             |     |      |
|--|----------------------------------|-----------------------------|-----|------|
| Hardcopy Color Image Storage                               | 1.2.840.10008.5.1.1.1.30         | See <a href="#">Table 5</a> | SCP | None |
| Hardcopy Grayscale Image Storage                           | 1.2.840.10008.5.1.1.1.29         | See <a href="#">Table 5</a> | SCP | None |
| Hemodynamic Waveform Storage                               | 1.2.840.10008.5.1.4.1.1.9.2.1    | See <a href="#">Table 5</a> | SCP | None |
| Key Object Selection Document                              | 1.2.840.10008.5.1.4.1.1.88.59    | See <a href="#">Table 5</a> | SCP | None |
| Mammography CAD SR   | 1.2.840.10008.5.1.4.1.1.88.50    | See <a href="#">Table 5</a> | SCP | None |
| Multi-frame Single Bit Secondary Capture Image Storage     | 1.2.840.10008.5.1.4.1.1.7.1      | See <a href="#">Table 5</a> | SCP | None |
| Multi-frame Grayscale Byte Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.2      | See <a href="#">Table 5</a> | SCP | None |
| Multi-frame Grayscale Word Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.3      | See <a href="#">Table 5</a> | SCP | None |
| Multi-frame True Color Secondary Capture Image Storage     | 1.2.840.10008.5.1.4.1.1.7.4      | See <a href="#">Table 5</a> | SCP | None |
| MR Image Storage   | 1.2.840.100008.5.1.4.1.1.4       | See <a href="#">Table 5</a> | SCP | None |
| Enhanced MR Image Storage                                  | 1.2.840.100008.5.1.4.1.1.4.1     | See <a href="#">Table 5</a> | SCP | None |
| Nuclear Medicine Image Storage                             | 1.2.840.10008.5.1.4.1.1.20       | See <a href="#">Table 5</a> | SCP | None |
| Nuclear Medicine Image Storage (Retired)                   | 1.2.840.10008.5.1.4.1.1.5        | See <a href="#">Table 5</a> | SCP | None |
| Ophthalmic Photography 8 Bit Image Storage                 | 1.2.840.10008.5.1.4.1.1.77.1.5.1 | See <a href="#">Table 5</a> | SCP | None |
| Ophthalmic Photography 16 Bit Image Storage                | 1.2.840.10008.5.1.4.1.1.77.1.5.2 | See <a href="#">Table 5</a> | SCP | None |
| Ophthalmic Tomography Image Storage                        | 1.2.840.10008.5.1.4.1.1.77.1.5.4 | See <a href="#">Table 5</a> | SCP | None |
| Positron Emission Tomography Image Storage                 | 1.2.840.10008.5.1.4.1.1.128      | See <a href="#">Table 5</a> | SCP | None |
| RT Beams Treatment Record Storage                          | 1.2.840.10008.5.1.4.1.1.481.4    | See <a href="#">Table 5</a> | SCP | None |
| RT Brachy Treatment Record Storage                         | 1.2.840.10008.5.1.4.1.1.481.6    | See <a href="#">Table 5</a> | SCP | None |
| RT Dose Storage  | 1.2.840.10008.5.1.4.1.1.481.2    | See <a href="#">Table 5</a> | SCP | None |
| RT Image Storage   | 1.2.840.10008.5.1.4.1.1.481.1    | See <a href="#">Table 5</a> | SCP | None |
| RT Plan Storage  | 1.2.840.10008.5.1.4.1.1.481.5    | See <a href="#">Table 5</a> | SCP | None |

|   |                                |                             |     |      |
|---|--------------------------------|-----------------------------|-----|------|
| RT Structure Set Storage                            | 1.2.840.10008.5.1.4.1.1.481.3  | See <a href="#">Table 5</a> | SCP | None |
| RT Treatment Summary Record Storage                 | 1.2.840.10008.5.1.4.1.1.481.7  | See <a href="#">Table 5</a> | SCP | None |
| Secondary Capture Image Storage                     | 1.2.840.10008.5.1.4.1.1.7      | See <a href="#">Table 5</a> | SCP | None |
| Stand-alone Curve Storage                           | 1.2.840.10008.5.1.4.1.1.9      | See <a href="#">Table 5</a> | SCP | None |
| Stand-alone Modality LUT Storage                    | 1.2.840.10008.5.1.4.1.1.10     | See <a href="#">Table 5</a> | SCP | None |
| Stand-alone Overlay Storage                         | 1.2.840.10008.5.1.4.1.1.8      | See <a href="#">Table 5</a> | SCP | None |
| Stand-alone VOI LUT Storage                         | 1.2.840.10008.5.1.4.1.1.11     | See <a href="#">Table 5</a> | SCP | None |
| Standalone PET Curve Storage                        | 1.2.840.10008.5.1.4.1.1.129    | See <a href="#">Table 5</a> | SCP | None |
| Stored Print Storage                                | 1.2.840.10008.5.1.1.27         | See <a href="#">Table 5</a> | SCP | None |
| Ultrasound Image Storage                            | 1.2.840.10008.5.1.4.1.1.6.1    | See <a href="#">Table 5</a> | SCP | None |
| Ultrasound Image Storage (Retired)                  | 1.2.840.10008.5.1.4.1.1.6      | See <a href="#">Table 5</a> | SCP | None |
| Ultrasound Multi-frame Image Storage                | 1.2.840.10008.5.1.4.1.1.3.1    | See <a href="#">Table 5</a> | SCP | None |
| Ultrasound Multi-frame Image Storage (Retired)      | 1.2.840.10008.5.1.4.1.1.3      | See <a href="#">Table 5</a> | SCP | None |
| VL Endoscopic Image Storage                         | 1.2.840.10008.5.1.4.1.1.77.1.1 | See <a href="#">Table 5</a> | SCP | None |
| VL Microscopic Image Storage                        | 1.2.840.10008.5.1.4.1.1.77.1.2 | See <a href="#">Table 5</a> | SCP | None |
| VL Slide-Coordinates Microscopic                    | 1.2.840.10008.5.1.4.1.1.77.1.3 | See <a href="#">Table 5</a> | SCP | None |
| VL Photographic Image Storage                       | 1.2.840.10008.5.1.4.1.1.77.1.4 | See <a href="#">Table 5</a> | SCP | None |
| X-Ray Angiographic Bi-Plane Image Storage (Retired) | 1.2.840.10008.5.1.4.1.1.12.3   | See <a href="#">Table 5</a> | SCP | None |
| X-Ray Angiographic Image Storage                    | 1.2.840.10008.5.1.4.1.1.12.1   | See <a href="#">Table 5</a> | SCP | None |
| X-Ray Radiofluoroscopic Image Storage               | 1.2.840.10008.5.1.4.1.1.12.2   | See <a href="#">Table 5</a> | SCP | None |

**Table 5. Transfer Syntax for Receive from a Remote System**

| Name  | UID                    |
|---|------------------------|
| Little Endian Explicit VR                               | 1.2.840.10008.1.2.1    |
| Little Endian Implicit VR                               | 1.2.840.10008.1.2      |
| Big Endian Explicit VR                                  | 1.2.840.10008.1.2.2    |
| Deflated Little Endian Explicit VR                      | 1.2.840.10008.1.2.1.99 |
| Lossless JPEG Image Compression (baseline)              | 1.2.840.10008.1.2.4.70 |
| Lossy JPEG Image Compression (8-bit, coding Process 1)  | 1.2.840.10008.1.2.4.50 |
| Lossy JPEG Image Compression (12-bit, coding Process 4) | 1.2.840.10008.1.2.4.51 |
| JPEG Lossless, Nonhierarchical (Processes 14)           | 1.2.840.10008.1.2.4.57 |
| JPEG 2000 Image Compression (Lossless Only)             | 1.2.840.10008.1.2.4.90 |
| JPEG 2000 Image Compression                             | 1.2.840.10008.1.2.4.91 |
| Lossless RLE Image Compression                          | 1.2.840.10008.1.2.5    |

## SOP Specific Conformance for all Storage SOP Classes

The iConnect Gateway DICOM Interface conforms to the SOPs of the Storage Service Class at Level 2 (full). No elements are discarded or coerced, unless explicitly configured to do so. Discarding offending DICOM elements (invalid tags) or coercing some data elements to address customer needs is possible; however, it can be done only with Merge Healthcare's direct involvement, and it implies custom implementations of coercion strategies.

In the event of a successful C-STORE operation, the object is fully stored to the iConnect Cloud Archive and is accessible through the iConnect Cloud Archive.

The iConnect Gateway performs minimal object validation to ensure database integrity. Additional validation may be configured to prevent the object from being stored in the archive.

Extended negotiation is not supported for Storage Service Classes.

iCloud Gateway returns Success status code (0x0000) when it successfully receives and saves object to disk for uploading to the iConnect Cloud Archive. Otherwise, failed response is returned with status code (0xC0110)

### *Presentation Context Acceptance Criterion*

There is no prioritization used for the acceptance of presentation contexts. Any combination

of supported transfer syntax and abstract syntax is accepted, in case the product is configured for them.

#### *Transfer Syntax Selection Policies*

The iConnect Gateway DICOM Interface selects transfer syntaxes in Acceptor first mode: the first transfer syntax in the archive's list of allowable transfer syntaxes that also exists in the proposed transfer syntax list will be selected.

## Response to a Query Request from a Remote System

#### *Associated Real-World Activity*

When the iConnect Gateway receives a query request (C-FIND) it will resolve the request against the iConnect Cloud Archive and return all found information to the query initiator. The gateway supports query transactions for all stored DICOM composite objects.

#### *Presentation Context Table*

The Enterprise DICOM Archive Manager will accept the following Transfer Syntaxes for each presentation context in an association it receives:

**Table 6. Presentation Context Table**

| Abstract Syntax                                       |                             | Transfer Syntax                    |                        | Role | Extended Negotiation |
|---|-----------------------------|------------------------------------|------------------------|------|----------------------|
| Name  | UID                         | Name                               | UID                    |      |                      |
| Patient Root Query/ Retrieve Information Model – FIND | 1.2.840.10008.5.1.4.1.2.1.1 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP  | Relational           |
|   |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |      |                      |
|   |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |      |                      |
|   |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |      |                      |
| Study Root Query/ Retrieve Information Model - FIND   | 1.2.840.10008.5.1.4.1.2.2.1 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP  | Relational           |
|   |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |      |                      |
|   |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |      |                      |
|   |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |      |                      |

|  |                             |                           |                        |     |            |
|--|-----------------------------|---------------------------|------------------------|-----|------------|
| Patient/<br>Study Only<br>Query/<br>Retrieve<br>Information<br>Model –<br>FIND | 1.2.840.10008.5.1.4.1.2.3.1 | Little Endian Explicit VR | 1.2.840.10008.1.2.1    | SCP | Relational |
|  |                             | Little Endian Implicit VR | 1.2.840.10008.1.2      |     |            |
|  |                             | Big Endian Explicit VR    | 1.2.840.10008.1.2.2    |     |            |
|  |                             | Deflated Little Endian    | 1.2.840.10008.1.2.1.99 |     |            |

#### *SOP Specific Conformance for all Find SOP Classes*

The iConnect Gateway conforms to the SOPs of the Query/Retrieve Service Class at both the hierarchical and relational level.

**Table 7. SOP Extended Negotiation**

| Item bytes | Field Name                   | Description of Field   |
|------------|------------------------------|--|
| 1          | Relational-queries/retrieval | This byte field defines relational-query/retrieval support for the Association-acceptor. It shall be encoded as an unsigned binary integer and shall use one of the following values<br>0 - relational-queries/retrieval not supported<br>1 - relational-queries/retrieval supported |

## Query Modes

The iConnect Gateway supports both the Hierarchical and Relational search modes.

## Information Models

The iConnect Gateway supports the Patient Root, Study Root, and Patient-Study Only information models.

**Table 8. Patient Root Image C-FIND Supported Attributes**

| Attribute Name       | Tag          | Query Keys Matching SCP | Query Keys Return SCP |
|----------------------|--------------|-------------------------|-----------------------|
| <b>Patient Level</b> |              |                         |                       |
| Patient Name         | (0010, 0010) | Supported               | Supported             |
| Patient ID           | (0010, 0020) | Supported (U)           | Supported             |

|                                     |              |               |           |
|-------------------------------------|--------------|---------------|-----------|
| Referenced Patient Sequence         | (0008, 1120) | Supported     | Supported |
| Referenced SOP Class UID            | (0008, 1150) | Supported     | Supported |
| Referenced SOP Instance UID         | (0008, 1155) | Supported     | Supported |
| Patient's Birth Date                | (0010, 0030) | Supported     | Supported |
| Patient's Birth Time                | (0010, 0032) | Supported     | Supported |
| Patient's Sex                       | (0010, 0040) | Supported     | Supported |
| Other Patient IDs                   | (0010, 1000) | Supported     | Supported |
| Other Patient Names                 | (0010, 1001) | Supported     | Supported |
| Ethnic Group                        | (0010, 2160) | Supported     | Supported |
| Patient Comments                    | (0010, 4000) | Supported     | Supported |
| Number of Patient Related Studies   | (0020, 1200) | Supported     | Supported |
| Number of Patient Related Series    | (0020, 1202) | Supported     | Supported |
| Number of Patient Related Instances | (0020, 1204) | Supported     | Supported |
| <b>Study Level</b>                  |              |               |           |
| Study Date                          | (0008, 0020) | Supported     | Supported |
| Study Time                          | (0008, 0030) | Supported     | Supported |
| Accession Number                    | (0008, 0050) | Supported     | Supported |
| Study ID                            | (0020, 0010) | Supported     | Supported |
| Study Instance UID                  | (0020, 000D) | Supported (U) | Supported |
| Modalities in Study                 | (0008, 0061) | Supported     | Supported |
| Referring Physician's Name          | (0008, 0090) | Supported     | Supported |
| Study Description                   | (0008, 1030) | Supported     | Supported |
| Procedure Code Sequence             | (0008, 1032) | Supported     | Supported |
| Code Value                          | (0008, 0100) | Supported     | Supported |
| Coding Scheme Designator            | (0008, 0102) | Supported     | Supported |
| Coding Scheme Version               | (0008, 0103) | Supported     | Supported |
| Code Meaning                        | (0008, 0104) | Supported     | Supported |
| Name of Physician(s) Reading Study  | (0008, 1060) | Supported     | Supported |
| Admitting Diagnoses Description     | (0008, 1080) | Supported     | Supported |

|                                     |              |               |           |
|-------------------------------------|--------------|---------------|-----------|
| Referred Study Sequence             | (0008, 1110) | Supported     | Supported |
| Referenced SOP Class UID            | (0008, 1150) | Supported     | Supported |
| Referenced SOP Instance UID         | (0008, 1155) | Supported     | Supported |
| Patient's Age                       | (0010, 1010) | Supported     | Supported |
| Patient's Size                      | (0010, 1020) | Supported     | Supported |
| Patient's Weight                    | (0010, 1030) | Supported     | Supported |
| Occupation                          | (0010, 2180) | Supported     | Supported |
| Additional Patient History          | (0010, 21B0) | Supported     | Supported |
| Other Study Numbers                 | (0020, 1070) | Supported     | Supported |
| Number of Study Related Series      | (0020, 1206) | Supported     | Supported |
| Number of Study Related Instances   | (0020, 1208) | Supported     | Supported |
| Interpretation Author               | (4008, 010C) | Supported     | Supported |
| <b>Series Level</b>                 |              |               |           |
| Modality                            | (0008, 0060) | Supported     | Supported |
| Series Number                       | (0020, 0011) | Supported     | Supported |
| Series Instance UID                 | (0020, 000E) | Supported (U) | Supported |
| Number of Series Related Instances  | (0020, 1209) | Supported     | Supported |
| Performed Procedure Step ID         | (0040, 0253) | Supported     | Supported |
| Reference Study Component Sequence  | (0008, 1111) | Supported     | Supported |
| Referenced SOP Class UID            | (0008, 1150) | Supported     | Supported |
| Referenced SOP Instance UID         | (0008, 1155) | Supported     | Supported |
| Request Attribute Sequence          | (0040, 0275) | Supported     | Supported |
| Requested Procedure ID              | (0040, 1001) | Supported     | Supported |
| Scheduled Procedure Step ID         | (0040, 0009) | Supported     | Supported |
| Performed Procedure Step Start Date | (0040, 0244) | Supported     | Supported |
| Performed Procedure Step Start Time | (0040, 0245) | Supported     | Supported |
| Body Part Examined                  | (0018, 0015) | Supported     | Supported |
| Station Name                        | (0008, 1010) | Supported     | Supported |

|  |              |               |           |
|--|--------------|---------------|-----------|
| Institution Name                       | (0008, 0080) | Supported     | Supported |
| Performing Physician's Name            | (0008, 1050) | Supported     | Supported |
| Manufacturer                           | (0008, 0070) | Supported     | Supported |
| Manufacturer's Model Name              | (0008, 1090) | Supported     | Supported |
| Series Description                     | (0008, 103E) | Supported     | Supported |
| Series Date                            | (0008, 0021) | Supported     | Supported |
| Series Time                            | (0008, 0031) | Supported     | Supported |
| Protocol Name                          | (0018, 1030) | Supported     | Supported |
| <b>Composite Object Instance Level</b> |              |               |           |
| Instance Number                        | (0020, 0013) | Supported     | Supported |
| Overlay Number                         | (0020, 0022) | Supported     | Supported |
| Curve Number                           | (0020, 0024) | Supported     | Supported |
| LUT Number                             | (0020, 0026) | Supported     | Supported |
| SOP Instance UID                       | (0008, 0018) | Supported (U) | Supported |
| SOP Class UID                          | (0008, 0016) | Supported     | Supported |
| <b>Image Specific Level</b>            |              |               |           |
| Rows                                   | (0020, 0010) | Supported     | Supported |
| Columns                                | (0020, 0011) | Supported     | Supported |
| Bits Allocated                         | (0028, 0100) | Supported     | Supported |
| Number of Frames                       | (0028, 0008) | Supported     | Supported |

**Table 9. Study Root Image C-FIND Supported Attributes**

| Attribute Name     | Tag          | Query Keys Matching SCP | Query Keys Return SCP |
|--------------------|--------------|-------------------------|-----------------------|
| <b>Study Level</b> |              |                         |                       |
| Study Date         | (0008, 0020) | Supported               | Supported             |
| Study Time         | (0008, 0030) | Supported               | Supported             |
| Accession Number   | (0008, 0050) | Supported               | Supported             |
| Patient Name       | (0010, 0010) | Supported               | Supported             |
| Patient ID         | (0010, 0020) | Supported (U)           | Supported             |



|                                    |              |               |           |
|------------------------------------|--------------|---------------|-----------|
| Study ID                           | (0020, 0010) | Supported     | Supported |
| Study Instance UID                 | (0020, 000D) | Supported (U) | Supported |
| Modalities in Study                | (0008, 0061) | Supported     | Supported |
| Referring Physician's Name         | (0008, 0090) | Supported     | Supported |
| Study Description                  | (0008, 1030) | Supported     | Supported |
| Procedure Code Sequence            | (0008, 1032) | Supported     | Supported |
| Code Value                         | (0008, 0100) | Supported     | Supported |
| Coding Scheme Designator           | (0008, 0102) | Supported     | Supported |
| Coding Scheme Version              | (0008, 0103) | Supported     | Supported |
| Code Meaning                       | (0008, 0104) | Supported     | Supported |
| Name of Physician(s) Reading Study | (0008, 1060) | Supported     | Supported |
| Admitting Diagnoses Description    | (0008, 1080) | Supported     | Supported |
| Referred Study Sequence            | (0008, 1110) | Supported     | Supported |
| Referenced SOP Class UID           | (0008, 1150) | Supported     | Supported |
| Referenced SOP Instance UID        | (0008, 1155) | Supported     | Supported |
| Referenced Patient Sequence        | (0008, 1120) | Supported     | Supported |
| Referenced SOP Class UID           | (0008, 1150) | Supported     | Supported |
| Referenced SOP Instance UID        | (0008, 1155) | Supported     | Supported |
| Referenced SOP Instance UID        | (0008, 1155) | Supported     | Supported |
| Patient's Birth Date               | (0010, 0030) | Supported     | Supported |
| Patient's Birth Time               | (0010, 0032) | Supported     | Supported |
| Patient's Sex                      | (0010, 0040) | Supported     | Supported |
| Other Patient IDs                  | (0010, 1000) | Supported     | Supported |
| Other Patient Names                | (0010, 1001) | Supported     | Supported |
| Patient's Age                      | (0010, 1010) | Supported     | Supported |
| Patient's Size                     | (0010, 1020) | Supported     | Supported |
| Patient's Weight                   | (0010, 1030) | Supported     | Supported |
| Ethnic Group                       | (0010, 2160) | Supported     | Supported |
| Occupation                         | (0010, 2180) | Supported     | Supported |

|                                     |              |               |           |
|-------------------------------------|--------------|---------------|-----------|
| Additional Patient History          | (0010, 21B0) | Supported     | Supported |
| Patient Comments                    | (0010, 4000) | Supported     | Supported |
| Other Study Numbers                 | (0020, 1070) | Supported     | Supported |
| Number of Patient Related Studies   | (0020, 1200) | Supported     | Supported |
| Number of Patient Related Series    | (0020, 1202) | Supported     | Supported |
| Number of Patient Related Instances | (0020, 1204) | Supported     | Supported |
| Number of Study Related Series      | (0020, 1206) | Supported     | Supported |
| Number of Study Related Instances   | (0020, 1208) | Supported     | Supported |
| Interpretation Author               | (4008, 010C) | Supported     | Supported |
| <b>Series Level</b>                 |              |               |           |
| Modality                            | (0008, 0060) | Supported     | Supported |
| Series Number                       | (0020, 0011) | Supported     | Supported |
| Series Instance UID                 | (0020, 000E) | Supported (U) | Supported |
| Number of Series Related Instances  | (0020, 1209) | Supported     | Supported |
| Performed Procedure Step ID         | (0040, 0253) | Supported     | Supported |
| Reference Study Component Sequence  | (0008, 1111) | Supported     | Supported |
| Referenced SOP Class UID            | (0008, 1150) | Supported     | Supported |
| Referenced SOP Instance UID         | (0008, 1155) | Supported     | Supported |
| Request Attribute Sequence          | (0040, 0275) | Supported     | Supported |
| Requested Procedure ID              | (0040, 1001) | Supported     | Supported |
| Scheduled Procedure Step ID         | (0040, 0009) | Supported     | Supported |
| Performed Procedure Step Start Date | (0040, 0244) | Supported     | Supported |
| Performed Procedure Step Start Time | (0040, 0245) | Supported     | Supported |
| Body Part Examined                  | (0018, 0015) | Supported     | Supported |
| Station Name                        | (0008, 1010) | Supported     | Supported |
| Institution Name                    | (0008, 0080) | Supported     | Supported |
| Performing Physician's Name         | (0008, 1050) | Supported     | Supported |
| Manufacturer                        | (0008, 0070) | Supported     | Supported |

|  |              |               |           |
|--|--------------|---------------|-----------|
| Manufacturer's Model Name              | (0008, 1090) | Supported     | Supported |
| Series Description                     | (0008, 103E) | Supported     | Supported |
| Series Date                            | (0008, 0021) | Supported     | Supported |
| Series Time                            | (0008, 0031) | Supported     | Supported |
| Protocol Name                          | (0018, 1030) | Supported     | Supported |
| <b>Composite Object Instance Level</b> |              |               |           |
| Instance Number                        | (0020, 0013) | Supported     | Supported |
| Overlay Number                         | (0020, 0022) | Supported     | Supported |
| Curve Number                           | (0020, 0024) | Supported     | Supported |
| LUT Number                             | (0020, 0026) | Supported     | Supported |
| SOP Instance UID                       | (0008, 0018) | Supported (U) | Supported |
| SOP Class UID                          | (0008, 0016) | Supported     | Supported |
| <b>Image Specific Level</b>            |              |               |           |
| Rows                                   | (0020, 0010) | Supported     | Supported |
| Columns                                | (0020, 0011) | Supported     | Supported |
| Bits Allocated                         | (0028, 0100) | Supported     | Supported |
| Number of Frames                       | (0028, 0008) | Supported     | Supported |

## Patient/Study Only Image C-FIND Supported Attributes

Attributes for the Patient and Study Levels of the Patient/Study Only Query/Retrieve Information Model are the same as the corresponding attributes for the Patient and Study Levels of the Patient Root Query/ Retrieve Information Model. For information about these attributes, see [Table 8](#).

**Table 10. Grayscale Soft Copy Presentation State C-FIND Supported Attributes**

| Attribute Name                      | Tag          | Query Keys Matching SCP | Query Keys Return SCP |
|-------------------------------------|--------------|-------------------------|-----------------------|
| <b>GSPS Instance Specific Level</b> |              |                         |                       |
| Presentation Label                  | (0070, 0080) | Supported               | Supported             |
| Presentation Description            | (0070, 0081) | Supported               | Supported             |

|                               |              |           |           |
|-------------------------------|--------------|-----------|-----------|
| Presentation Creation Date    | (0070, 0082) | Supported | Supported |
| Presentation Creation Time    | (0070, 0083) | Supported | Supported |
| Presentation Creator's Name   | (0070, 0084) | Supported | Supported |
| Referenced Series Sequence    | (0008, 1115) | Supported | Supported |
| >Series Instance UID          | (0020, 000E) | Supported | Supported |
| >Referenced Image Sequence    | (0008, 1140) | Supported | Supported |
| >>Referenced SOP Class UID    | (0008, 1150) | Supported | Supported |
| >>Referenced SOP Instance UID | (0008, 1155) | Supported | Supported |

**Table 11. Table 15: Basic Text SR and Enhanced SR C-FIND Supported Attributes**

| Attribute Name                                   | Tag          | Query Keys Matching SCP | Query Keys Return SCP |
|--|--------------|-------------------------|-----------------------|
| <b>SR Instance Specific Level</b>                |              |                         |                       |
| Completion Flag                                  | (0040, A491) | Supported               | Supported             |
| Verification Flag                                | (0040, A493) | Supported               | Supported             |
| Content Date                                     | (0008, 0023) | Supported               | Supported             |
| Content Time                                     | (0008, 0033) | Supported               | Supported             |
| Observation Date Time                            | (0040, A032) | Supported               | Supported             |
| Verifying Observer Sequence                      | (0040, A073) | Supported               | Supported             |
| >Verifying Organization                          | (0040, A027) | Supported               | Supported             |
| >Verification DateTime                           | (0040, A030) | Supported               | Supported             |
| >Verifying Observer Name                         | (0040, A075) | Supported               | Supported             |
| >Verifying Observer Identification Code Sequence | (0040, A088) | Supported               | Supported             |
| Referenced Request Sequence                      | (0040, A370) | Supported               | Supported             |
| >Study Instance UID                              | (0020, 000D) | Supported               | Supported             |
| >Accession Number                                | (0008, 0050) | Supported               | Supported             |
| >Requested Procedure ID                          | (0040, 1000) | Supported               | Supported             |
| >Requested Procedure Code Sequence               | (0032, 1064) | Supported               | Supported             |
| >>Code Value                                     | (0008, 0100) | Supported               | Supported             |

|                            |              |           |           |
|----------------------------|--------------|-----------|-----------|
| >>Code Scheme Designator   | (0008, 0102) | Supported | Supported |
| >>Code Scheme Version      | (0008, 0103) | Supported | Supported |
| >>Code Meaning             | (0008, 0104) | Supported | Supported |
| Concept Name Code Sequence | (0040, A043) | Supported | Supported |
| >Code Value                | (0008, 0100) | Supported | Supported |
| >Coding Scheme Designator  | (0008, 0102) | Supported | Supported |
| >Coding Scheme Version     | (0008, 0103) | Supported | Supported |
| >Code Meaning              | (0008, 0104) | Supported | Supported |

---

## C-FIND Status Codes

---

If the Enterprise DICOM Archive Manager returns one of the following status codes for C-FIND then the operation was unsuccessful. All status codes are in hexadecimal.

C001 (Unable To Process) - Indicates that the archive was unable to fully process the query. The Error Comment field of the status will be populated with a description of the error encountered.

---

## Query Implementation Specifics

---

The SCUs are recommended to append wildcard "\*" at the end of each component of any structured name to facilitate matching.

### *Presentation Context Acceptance Criterion*

There is no prioritization used for the acceptance of presentation contexts. Any combination of supported transfer syntax and abstract syntax is accepted, in case the product is configured for them.

### *Transfer Syntax Selection Policies*

The iConnect Gateway selects transfer syntaxes in Acceptor first modes: the first transfer syntax in the archive's list of allowable transfer syntaxes that also exists in the proposed transfer syntax list will be selected.

---

## Respond to a Retrieve Request from a Remote System

---

### *Associated Real-World Activity*

When the iConnect Gateway receives a request for retrieval (C-MOVE), the request will be resolved against the iConnect Cloud Archive's underlying database. If any objects are found, it will attempt to establish an association over which to store the found objects (C-MOVE only). If the association is established, it will retrieve the found objects from the underlying platform and store them over the association, returning pending responses to the retrieval request initiator.

#### Presentation Context Table

The iConnect Gateway will accept the following Transfer Syntaxes for each presentation context in an association it receives:

**Table 12. Presentation Context Table**

| Abstract Syntax  |                             |                                    |                        | Role | Extended Negotiation |
|--|-----------------------------|------------------------------------|------------------------|------|----------------------|
| Name   | UID                         | Name                               | UID                    |      |                      |
| Patient Root Query/Retrieve Information Model – MOVE       | 1.2.840.10008.5.1.4.1.2.1.2 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP  | Relational           |
|  |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |      |                      |
|  |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |      |                      |
|  |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |      |                      |
| Patient Root Query/Retrieve Information Model – GET        | 1.2.840.10008.5.1.4.1.2.1.3 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP  | Relational           |
|  |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |      |                      |
|  |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |      |                      |
|  |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |      |                      |
| Study Root Query/Retrieve Information Model – MOVE         | 1.2.840.10008.5.1.4.1.2.2.2 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP  | Relational           |
|  |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |      |                      |
|  |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |      |                      |
|  |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |      |                      |
| Study Root Query/Retrieve Information Model – GET          | 1.2.840.10008.5.1.4.1.2.2.3 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP  | Relational           |
|  |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |      |                      |
|  |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |      |                      |
|  |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |      |                      |
| Patient/Study Only Query/Retrieve Information Model – MOVE | 1.2.840.10008.5.1.4.1.2.3.2 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP  | Relational           |
|  |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |      |                      |
|  |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |      |                      |
|  |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |      |                      |

|   |                             |                                    |                        |     |            |
|---|-----------------------------|------------------------------------|------------------------|-----|------------|
| Patient/Study Only Query/Retrieve Information Model - GET | 1.2.840.10008.5.1.4.1.2.3.3 | Little Endian Explicit VR          | 1.2.840.10008.1.2.1    | SCP | Relational |
|   |                             | Little Endian Implicit VR          | 1.2.840.10008.1.2      |     |            |
|   |                             | Big Endian Explicit VR             | 1.2.840.10008.1.2.2    |     |            |
|   |                             | Deflated Little Endian Explicit VR | 1.2.840.10008.1.2.1.99 |     |            |

*SOP Specific Conformance for All Query/Retrieve SOP Classes*

## C-MOVE Status Codes

If the iConnect Gateway returns one of the following status codes for C-MOVE then the operation was unsuccessful. All status codes are in hexadecimal.

C001 (Unable To Process) - Indicates that the archive was unable to fully process the move/get. The Error

Comment field of the status will be populated with a description of the error encountered.

A801 (Move Destination Unknown) - Indicates that the archive was unable to connect to the move destination (C-MOVE only).

### *Presentation Context Acceptance Criterion*

There is no prioritization used for the acceptance of presentation contexts. Any combination of supported transfer syntax and abstract syntax is accepted, in case the product is configured for them.

### *Transfer Syntax Selection Policies*

The iConnect Gateway selects transfer syntaxes in Acceptor first mode: the first transfer syntax in the archive's list of allowable transfer syntaxes that also exists in the proposed transfer syntax list will be selected.

# Configuration

The iConnect Gateway tool (a graphical user interface (GUI)) is used for configuration and administration. This tool is remotely enabled through the use of Java Remote Method Invocation (RMI), allowing the system to be configured and administrated from any location on the network. Configurable parameters include, but are not limited to, system level debugging, known remote application entities, various security concerns, and service classes supported.

## AE Title/Presentation Address Mapping

All AE Title/Presentation Address mapping may need to be present in multiple locations, depending on the Application Entity. In all cases, iConnect Gateway allows for the configuration of the AE Title.

## Configurable Parameters

The amount of configuration the system allows is extensive.

In general, the following parameters are configurable for each component:

- Listening IP port number
- Application entity title
- Number of simultaneous associations
- DICOM Part 10 file storage location
- Accepted SOP Classes
- Remote AE list (AE Title, Hostname, Port)
- Network Timeout

A member of Merge Healthcare's support staff must perform configuration.