

Merge Universal Viewer

9.2

IHE INTEGRATION STATEMENT

© Copyright Merge Healthcare Solutions Inc. 2026.

Licensed materials – Property of Merge Healthcare Solutions Inc.

The content of this document is confidential information of Merge Healthcare Solutions Inc. and its use and disclosure is subject to the terms of the agreement pursuant to which you obtained the software that accompanies the documentation.

Merge Healthcare and the Merge Healthcare logo are trademarks of Merge Healthcare Inc.

Microsoft, Windows, and the Windows logo are trademarks of Microsoft Corporation in the United States, other countries, or both. All other names are trademarks or registered trademarks of their respective companies.

U.S. GOVERNMENT RESTRICTED RIGHTS:

This product is a “Commercial Item” offered with “Restricted rights.” The Government’s right to use, modify, reproduce, release, perform, display or disclose this documentation are subject to the restrictions set forth in Federal Acquisition Regulation (“FAR”) 12.211 and 12.212 for civilian agencies and in DFARS 227.7202-3 for military agencies. Contractor is Merge Healthcare Solutions Inc.

Manufacturer:

Merge Healthcare Incorporated, an affiliated of Merative US, L.P
900 Walnut Ridge Drive
Hartland, WI 53029 USA

For application support or to report issues with user documentation, contact Customer Support:

☎ 1-877-741-5369 (North America)
+44203808.3608 (Europe, the Middle East and Africa)

✉ mergesupport@merative.com

Part	Date	Revision	Description
MUV-SMD-0001	July 2026	01	Initial revision for Merge Universal Viewer 9.2

IHE Integration Statement		Date
		2025
Vendor	Product Name	Version
Merge Healthcare Solutions Inc.	Merge Universal Viewer	9.2
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Audit Trail and Node Authentication	Secure Application	Radiology Audit trail, FQDN Validation of Server Certificate
Consistent Time	Time Client	None
Patient Demographics Query	Patient Demographics Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Cross-Community Access	Document Consumer	None
Cross-Community Access for Imaging	Imaging Document Consumer	None
Cross-Enterprise Document Sharing	Document Consumer	None
Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer	None
Cross-Enterprise Document Sharing of Scanned Documents	Content Consumer	View Option
Internet address for vendor’s IHE information: https://www.merative.com/support/merge/merge-healthcare-statements		
Links to Standard Conformance Statements for the Implementation		
HL7	None	
DICOM	https://www.merative.com/support/merge/merge-healthcare-statements	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org