



IBM iConnect[®] Enterprise Archive

14.2

IHE Integration Statement

Licensed Materials - Property of IBM

© Copyright IBM Corporation 2022.

IBM, the IBM logo, IBM Watson, and Watson Health are trademarks of IBM Corporation in the United States, other countries or both. The use of IBM trademarks herein are by authorization and license from the IBM Corporation.

Microsoft, Windows, Windows NT, and the Windows logo are trademarks of Microsoft Corporation in the United States, other countries, or both.

All other company or product names are trademarks of their respective companies.

Licensees of this program who wish to have information about it for the purpose of enabling: (i) the exchange of information between independently created programs and other programs (including this one) and (ii) the mutual use of the information which has been exchanged, should contact the IBM Director of Licensing. Such information may be available, subject to appropriate terms and conditions, including in some cases, payment of a fee.

IBM Director of Licensing

IBM Corporation

North Castle Drive, MD-NC119

Armonk, NY 10504-1785

IBM grants limited permission to licensees to make hardcopy or other reproductions of any machine-readable documentation, provided that each such reproduction shall carry the IBM copyright notices and that use of the reproduction shall be governed by the terms and conditions specified by IBM in the Licensed Program Specifications. Any reproduction or use beyond the limited permission granted herein shall be a breach of the license agreement and an infringement of the applicable copyrights.

Materials displayed or reproduced by this program may be protected by copyright and/or contract restrictions of IBM and/or others. The user is responsible for having permission to display and reproduce such materials and for including applicable copyright notices

For application support or to report issues with user documentation, contact Customer Support:

☎ 1-877-741-5369 (North America)
+31.20.514.5073 (Europe, the Middle East, and Africa)
1800 316 746 (Australia)

✉ WHISupport@us.ibm.com

Part	Date	Revision	Description
ECM-41366	May 2022	1.0	Initial version for 14.2

The latest version of this document can be found at <https://mergecustomer.force.com/mergeusercommunity/login>.

Chapter 1. IHE Integration Statement

The following table outlines the IHE Integration Statement.

IHE Integration Statement		Date
		2022
Vendor	Product Name	Version
Merge Healthcare Incorporated	IBM iConnect® Enterprise Archive	14.2
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Patient Information Reconciliation	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Consistent Presentation of Images	Image Manager / Image Archive	None
Presentation of Grouped Procedures	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Evidence Documents	Image Manager / Image Archive	None
Key Image Note	Image Manager / Image Archive	None
Simple Image and Numeric Report	Report Repository	None
	Enterprise Report Repository	None
	External Report Repository Access	None
Access to Radiology Information	Image Manager / Image Archive	None
	Report Repository	None
	External Report Repository Access	None
Nuclear Medicine	Image Manager / Image Archive	None
Cardiac Catheterization Workflow	Image Manager / Image Archive	Intermittently Connected Modality, Cardiac Cath
	Performed Procedure Step Manager	None
Echo Workflow	Image Manager / Image Archive	Echocardiography, Intermittently Connected Modality
	Performed Procedure Step Manager	None
Evidence Documents - Cardiology	Image Manager / Image Archive	None
Consistent Time (CT)	Time Client	None

Mammography	Image Manager / Image Archive	None
Stress Testing Workflow	Image Manager / Image Archive	None
	Performed Procedure Step Manager	None
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Image Document Source	Set of DICOM Instances
	Document Repository	None
Cross-Enterprise Document Sharing (XDS.b)	Document Registry	Patient Identity Feed (HL7 V2), Patient Identity Feed (HL7 V3), On-Demand Documents, Reference ID, Document Metadata Update
	Document Repository	None
	Integrated Document Source/Repository	Document Replace
Audit Trail and Node Authentication (ATNA)	Secure Application	None
Patient Identifier Cross-reference (PIX)	Patient Identity Cross-reference Manager	None
	Patient Identifier Cross-reference Consumer	PIX Update Notification
Patient Demographics Query (PDQ)	Patient Demographics Supplier	None
	Patient Demographics Consumer	None
Imaging Object Change Management (IOCM)	Image Manager/Archive	None
	Change Requester	None
Document Metadata Subscription (DSUB)	Document Metadata Notification Broker	Patient-Independent Subscription, Folder Subscription
	Document Metadata Publisher	None
Multi-Patient Queries (MPQ)	Document Registry	None
Web-based Image Capture (WIC)	Image Manager	None
Internet address for vendor's IHE information:	www.merge.com/Support/Resources/IHE-Statements.aspx	
Links to Standards Conformance Statements for the Implementation		
HL7	Refer to Merge Communities -- http://www.merge.com/Support/Support-Portal.aspx	
DICOM	www.merge.com/Support/Resources/DICOM-Conformance-Statements.aspx	
Links to general information on IHE		
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.net	In Japan: http://www.ihe-j.org