




V11.2 UPDATES

Merge Hemo

Merge Hemo

Merge Hemo™, a KLAS Category Leader for Hemodynamics for over a decade, and Best in KLAS for 2025, simplifies data collection by mimicking your existing workflow. The intuitive touch screen and flexible interface facilitate the rapid documentation of study data, making your cath lab more streamlined, efficient, and cost-effective.

Merge delivers a comprehensive and intuitive solution to support your entire cath lab allowing for more accurate and streamlined data capture.



Your flexible transition path to the future of hemodynamics

Merge 11.2 offers a mixed mode for customers to upgrade from their 10.x hardware and software to the new 11.2 system over time.

Merge values your investment in our products, and we want you to upgrade on your own terms, in a way that focuses on your ROI. Inquire with your Merge partner about creating a hybrid environment for your Merge Hemo station.

What's new for Hemo 11.2

Upgrading to 11.2 is essential for organizations seeking greater reliability and operational efficiency in hemodynamic monitoring. Included in this release is:

- **Windows 11 LTSC Compatibility:** Direct support for Windows 11 LTSC ensures a future-proof infrastructure, enhanced security, and compliance with IT standards.
- **Nellcor SpO2 Compatibility:** Adds support for Nellcor SpO2 monitoring devices, enabling more flexible, comprehensive patient monitoring and streamlined workflows.
- **Streamlined DC Workstation Deployment:** Allows IT to push mass updates to all DC workstations from the server, reducing update time across multiple labs and facilities.
- **Updated Technical Requirements:** Includes compatibility with Microsoft Windows Server 2019/2022, SQL Server 2022, and prerequisite software like OpenJDK 23.0.1 and MirthConnect 4.5.0.

With your Merge Hemo upgrade you will experience:

- **Enhanced Security:** Strengthens system resilience through proactive enhancements, ensuring compliance with current IT standards and hardware requirements.
- **Support for New Operating Systems:** Addresses OS hardware end-of-life risks by supporting the latest operating systems and hardware.
- **Improved System Stability:** Remediates known issues from previous releases, enhancing reliability for clinical workflows.

About Merge

Trusted by 6 of the 10 largest U.S. health systems, Merge by Merative empowers healthcare organizations with advanced medical imaging solutions to enhance workflows, optimize care delivery, and help improve patient outcomes. Merge's Imaging Suite, built on a cloud-native foundation, provides solutions for Vendor-Neutral Archive (VNA), PACS, Enterprise Image Viewing and Workflow Orchestration. Merge's portfolio also includes Best in KLAS Cardiology and Hemodynamic Monitoring, and Digital Pathology.

Learn more at merative.com/merge-imaging

About Merative

Merative provides data, analytics, and software for healthcare and government social services. With focused innovation and deep expertise, Merative works with providers, employers, health plans, governments, and life sciences companies to drive real progress. Merative helps clients orient information and insights around the people they serve to improve decision-making and performance.

Learn more at merative.com

© Merative US L.P. 2025. All Rights Reserved.

Produced in the United States of America
November 2025

Merative and the Merative logo are trademarks of Merative US L.P. Other product and service names might be trademarks of Merative or other companies.

The information contained in this publication is provided for informational purposes only. While efforts were made to verify the completeness and accuracy of the information contained in this publication, it is provided AS IS without warranty of any kind, express or implied. In addition, this information is based on Merative's product plans and strategy as of the date of this publication, which are subject to change by Merative without notice. Nothing contained in this publication is intended to, nor shall have the effect of, creating any warranties or representations from Merative, or stating or implying that any activities undertaken by you will result in any specific performance results. Merative products are warranted according to the terms and conditions of the agreements under which they are provided.

Merge imaging solutions are manufactured by Merge Healthcare Incorporated, an affiliate of Merative US, L.P.
900 Walnut Ridge Drive, Hartland, WI USA 53029.

HEMO-44770 Rev 1.0

Learn how you can benefit

Merative provides award-winning enterprise imaging solutions for radiology, cardiology and orthopedics, using vendor neutral cloud and AI technologies to reduce workspace complexity and create seamless workflows.

Visit [our website](#) or [talk to an enterprise imaging solutions expert](#) to learn more about how Merative can help you.

