



Merge Cardio™

12.2

IHE INTEGRATION STATEMENT

**© Copyright Merge Healthcare Solutions Inc. 2022.**

The content of this document is confidential information of Merge Healthcare Solutions Inc. and its use and disclosure is subject to the terms of the agreement pursuant to which you obtained the software that accompanies the documentation.

Merge Healthcare and the Merge Healthcare logo are trademarks of Merge Healthcare Inc.

Microsoft, Windows, and the Windows logo are trademarks of Microsoft Corporation in the United States, other countries, or both.

All other names are trademarks or registered trademarks of their respective companies.

**U.S. GOVERNMENT RESTRICTED RIGHTS:**

This product is a “Commercial Item” offered with “Restricted Rights.” The Government's rights to use, modify, reproduce, release, perform, display or disclose this documentation are subject to the restrictions set forth in Federal Acquisition Regulation (“FAR”) 12.211 and 12.212 for civilian agencies and in DFARS 227.7202-3 for military agencies. Contractor is Merge Healthcare Solutions Inc.

For application support or to report issues with user documentation, contact Customer Support:

☎ 1-877-741-5369 (North America)

✉ [mergesupport@merative.com](mailto:mergesupport@merative.com)

Part	Date	Revision	Description
VRCS-35508	June 2022	1.0	Release 12.2.

The latest version of this document can be found at <https://mergecustomer.force.com/mergeusercommunity/login>.

<b>IHE Integration Statement</b>		<b>Date</b>
		2022
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
Merge Healthcare Solutions Inc.	Merge Cardio	12.2
This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors, and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
<b>Scheduled Workflow</b>	Image Manager / Image Archive	
	Performed Procedure Step Manager	
	Image Display	
	Evidence Creator	Creator Performed Procedure Step
<b>Patient Information Reconciliation</b>	Image Manager / Image Archive	
	Performed Procedure Step Manager	
<b>Consistent Presentation of Images</b>	Image Manager / Image Archive	
	Image Display	
	Evidence Creator	
	Print Composer	
<b>Presentation of Grouped Procedures</b>	Image Manager / Image Archive	
	Performed Procedure Step Manager	
<b>Evidence Documents</b>	Image Manager / Image Archive	
	Image Display	
<b>Key Image Note</b>	Image Manager / Image Archive	
	Image Display	
	Evidence Creator	
<b>Simple Image and Numeric Report</b>	Report Repository	
	Report Reader	
	Enterprise Report Repository	
	External Report Repository Access	
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
<b>Access to Radiology Information</b>	Image Manager / Image Archive	Multiple Sources
	Report Repository	Multiple Sources
	Image Display	Multiple Sources
	Report Reader	Multiple Sources
	External Report Repository Access	

<b>Nuclear Medicine</b>	Image Manager / Image Archive	
<b>Cardiac Catheterization Workflow</b>	Image Manager / Image Archive	Intermittently Connected Modality Cardiac Cath
	Performed Procedure Step Manager	
<b>Echo Workflow</b>	Image Manager / Image Archive	Echocardiography Intermittently Connected Modality
	Performed Procedure Step Manager	
<b>Evidence Documents - Cardiology</b>	Image Manager / Image Archive	
<b>Retrieve ECG for Display</b>	Display	
<b>Retrieve Information for Display</b>	Display	
<b>Consistent Time</b>	Time Client	
<b>Patient Synchronized Applications</b>	Patient Context Participant	
<b>Portable Data for Imaging</b>	Portable Media Creator	
	Portable Media Importer	
<b>Mammography</b>	Image Manager / Image Archive	
<b>Stress Testing Workflow</b>	Image Manager / Image Archive	
	Performed Procedure Step Manager	
<b>Internet address for vendor's IHE information:</b> <a href="https://www.ibm.com/support/pages/node/6262447#IHE">https://www.ibm.com/support/pages/node/6262447#IHE</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b>	<a href="https://www.ibm.com/support/pages/node/6262447#DICOM">https://www.ibm.com/support/pages/node/6262447#DICOM</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>